fluorocare® SARS-CoV-2 Antigen Test Kit

CE 2934

PRODUCT NAME

Common Name: **SARS-CoV-2 Antigen Test Kit(Colloidal Gold Chromatographic Immunoassay) (Saliva)**

Product Code: **MF-91**

PACKING

1 Test/box, 2 Tests/box, 5 Tests/box

INTENDED USE

The fluorocare® SARS-CoV-2 Antigen Test Kit is applicable to the qualitative detection novel Coronavirus (SARS-CoV-2) Antigen in population saliva samples in vitro.

It is intended as an aid in the diagnosis of coronavirus infection disease (COVID-19) for symptomatic patients within 7 days after onset of symptoms, which is caused by SARS-CoV-2.

For in vitro diagnostic use only. For self-testing use. This kit is suitable for people over 2 years old. People under the age of 2-12 cannot operate by themselves. This kit should be used by adults or parents (18-60 years old) for sample collection and testing.

People aged 13-17 can use this kit to collect samples and test samples under the supervision of adults or parents (18-60 years old). Supervisors should ensure that users have a detailed understanding of the requirements of the instructions and watch whether the user's operation is correct.

For people over 75 years old, it is recommended that family members or guardians (18-60 years old) use this kit to collect samples and test samples.

INTRODUCTION

The novel coronaviruses belong to the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection. Asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

PRINCIPLE

The novel Coronavirus (SARS-CoV-2) Antigen test is a qualitatively test to detect SARS-CoV-2 Antigen in saliva samples by the colloidal gold method. After pre-treat the samples, the novel Coronavirus (SARS-CoV-2) Antigen in the sample to be tested is combined with the novel coronavirus (SARS-CoV-2) antibody labeled with colloidal gold on the Conjuagte pad to form SARS-CoV-2 Antigen-SARS-CoV-2 antibody-colloidal gold complex. Due to chromatography, the SARS-CoV-2 Antigen-antibody complex binds to the antibody enclosed within the detection line area, the SARS-CoV-2 Antigen-antibody complex binds to the antibody enclosed within the detection line area, showing a purple band. Colloidal gold-labelled SARS-CoV-2 antibody diffused to the quality control line (C) region and is captured by Goat anti-mouse IgG to form purple bands. When the reaction is over, the results can be judged by visual observation.

Specifications

Components	Specifications			
	REF F-91-1 1 Test/box	REF MF-91-2 2 Tests/box	REF F-91-5 5 Tests/box	Components
Test Card (including the desiccant)	1 cassette	2 cassette	5 cassettes	Each test card is mainly composed of a plastic casing and a test strip. The test strip contains: The Nitrocellulose membrane is coated with SARS-CoV-2 antibody, and the Conjugate pad contains colloidal gold-labelled SARS-CoV-2 antibody, and other components include PVC pad and absorbent paper.
Instruction of use	1 copy	1 copy	1 copy	/
Saliva collection cup	1 piece	2 piece	5 pieces	/
Prefilled Sample treatment tube	1 piece	2 piece	5 pieces	The bottle contains buffers of PBS, Trolaton-100, NP-40, and SDS. 1 mL per tube.

WHAT ELSE DO YOU NEED?

Timer.

STORAGE CONDITION AND EXPIRY DATE

1. Test kit store at 2-30°C in dry place and protect from light. Test kit is valid for 18 months.

2. The Test Card must remain in the sealed pouch until use.




3. The test reagent card is valid for use within 60 minutes after opening the package

HOW TO USE THE TEST?


Use a disinfectant to disinfect your hands after washing your hands, and the test needs to be done on a clean table. Before testing, read the operating instructions carefully, and restore the testing kit and samples to room temperature (20- 25°C) before using. The test should be done at 20-25°C. If the kit is removed from the refrigerator, allow it to stand at room temperature(20- 25°C) for 5 minutes before testing.

1. Sample collection

Saliva collection method:

4. Install the dropper lid of sample treatment tube and shake the tube up and down for 5 times. Please be careful to avoid liquid leakage. Then wait for the reaction between the saliva sample and the treatment solution for 1 minute before proceeding to the next stage of the test.

5. After 1 minute of sample reaction time, shake the tube up and down 5 times ,Please be careful to avoid liquid leakage. Unscrew the terminal at the top of the cap. If the terminal on the top of the cap is not unscrewed, and if the blue stopper inside the sample processing tube is not removed, it will not be possible to drip liquid!

And vertically add 2 drops (approximately 60 μ L) of sample solution to the sample hole(Marked as "S").Only 2 drops of the treated sample solution can be added ! Adding too much or too little of the treated sample solution may result in invalid test results !

Ensure the test remains on a flat surface and at room temperature for 15 minutes.

After the sample has been added, the cap, the top terminal of the cap and the blue stopper are all capped back into the treatment tube and treated as contaminants.

home to avoid spreading the virus to others, wear a face mask when recommended and wash your hands regularly with soap and water.

2. A negative result means the virus that causes COVID-19 was not found in your sample. A negative test result does not guarantee that you do not or have never had COVID-19, nor does it confirm whether or not you are currently contagious. If you have cold symptoms, dyspnea or high fever, you should assume that you have covid-19 because the home test does not provide complete certainty

You can contact your doctor to find out if another test is needed. In the meantime, try to avoid leaving your home and have as little contact as possible with others, including the people you live with. Use disposable tissues and throw them straight in the bin. Sneeze and cough into the crook of your elbow. Wash your hands regularly and wear a face mask.

3. You should not make any medically related decisions without first consulting your doctor. Actions you take after getting your test results must comply with current local regulations.

LIMITATION OF METHODOLOGY

1. This kit is a qualitative test and is only used for in vitro auxiliary diagnosis.

2. Negative test results may occur if the level of antigen in a sample is below the detection limit of the test, or from improper sample collection, and the negative results are not intended to exclude other non 2019-nCoV virus infections.

3. Unreasonable sampling, transportation, handling, and low virus content in samples may lead to false negatives.

4. This reagent is a qualitative assay. As it is with any diagnostic procedure, a confirmed 2019-nCoV infection diagnosis should only be made by a physician after evaluating all clinical and laboratory findings.

5. Reading the test results earlier than 15 minutes or later than 20 minutes may give incorrect results.

6. A negative test result does not rule out a coronavirus infection and does not exempt you from the applicable rules for spread control (e.g. contact restrictions and protective measures).

QUESTION & ANSWER

Q1. How does the SARS-CoV-2 Antigen Test Kit (Colloidal Gold Chromatographic Immunoassay) work?

The SARS-CoV-2 Antigen Test Kit is an antigen test that is to detect novel Coronavirus (SARS-CoV-2) Antigen in population saliva samples in vitro.

Q2.What is the difference between a COVID-19 antigen, molecular, and antibody test?

There are different kinds of tests for diagnosing COVID-19. Molecular tests (also known as PCR tests) detect genetic material from the virus.

Antigen tests are very specific for the virus but are not as sensitive as molecular tests.

Another type of test is an antibody test. A COVID-19 antibody test detects antibodies that have been made by your immune system in response to a previous COVID-19 infection.

Q3. Will this test hurt?

No, Just collect saliva, do not cause harm.

Q4. What does it mean if I have a positive test result?

A positive result means that you may have COVID-19 disease. Please contact your doctor for further medical suggestion. It is likely you will be asked to isolate yourself at home to avoid spreading the virus to others, wear a face mask when recommended and wash your hands regularly with soap and water. A positive result does not in any way guarantee that you are or will be immune and therefore cannot (or can no longer) become infected.


You should not make any medically related decisions without first consulting your doctor.

If there is a infection of COVID-19 virus,You should pay attention to personal protection to prevent infecting others, and go to the hospital for diagnosis as soon as possible.

Q5. What does it mean if I have a negative test result?

A negative result means the virus that causes COVID-19 was not found in your sample.

A negative test result does not guarantee that you do not or have never had COVID-19, nor does it confirm whether or not you are currently contagious. Do you have cold symptoms in addition to the negative at-home test? Since the at-home test does not provide complete certainty, you should assume that you have COVID-19. You can contact your doctor to find out if another test is needed. In the meantime, try to avoid leaving your home and have as little contact as possible with others, including the people you live with. Use disposable tissues and throw them straight in the bin. Sneeze and cough into the crook of your elbow. Wash your hands regularly and wear a face mask. Are your symptoms getting worse (difficulty breathing, high fever, etc.)? Contact your doctor/health


6. The test card is kept at room temperature for 15 minutes to observe the test results, but the observation results over 20 minutes were invalid.So read the results within 15-20 minutes after dropping the sample. Do not touch or lift the card from a flat surface during the course of the test.

TAKE MEASURES DURING THE TEST TO PREVENT SPREAD INFECTION

1. After the completion of observation and testing, put the used product components into a plastic bag, close and put the bag into another plastic bag and discard it. Reapply hand sanitizer to disinfect your hands.

2. Please complete the above test operation alone in an isolated room.


HOW TO READ THE RESULTS?

1. Positive: Two purple lines, both the detection line (T line) and the quality control line (C line) display color.

NOTE: It does not matter the line (T) is lighter or darker than the other; the result is "Positive".

2. Negative: a purple line, quality control line (C line) display color.

3. Invalid: The position of the quality control line (C line) in the observation window does not show any color rendering. These indicate that the test is invalid, so the sample should be re-sampled for testing.

WHAT SHOULD YOU DO AFTER READING THE TEST RESULT ?

1. A positive result means that you may have COVID-19 disease. Please contact your doctor for further medical suggestion. It is likely you will be asked to isolate yourself at

provider immediately.

You should not make any medically related decisions without first consulting your doctor.

Actions you take after getting your test results must comply with current local regulations.

Q6. How accurate is the SARS-CoV-2 Antigen Test Kit (Colloidal Gold Chromatographic Immunoassay)?

The test has been shown in field clinical evaluations performed by professional health care persons to correctly identify 96.86% (792/818) of 2019-nCoV samples (known as the test's accuracy). Further, in field clinical evaluations conducted, the test correctly identified 100% (100/100) 2019-nCoV negative samples when performed by self test users.

Q7. Is there any chance that I get a "false" negative result with this test?

It is possible for this test to give an incorrect negative (false negative) result". This means that you could still have COVID-19 even though the test result is negative. If your result is negative and you still experience symptoms related to COVID-19 , such as fever, cough and/or shortness of breath, you should seek help from your healthcare provider.

Q8. Is there any chance that I get an incorrect positive result?

There is a very small chance that this test gives you a positive result that is incorrect (false positive). If you get a positive result, you should self-isolate and seek medical help from you healthcare provider.

Q9. I have used the test but no colored band appears at control line (C). What should I do?

If there is no colored band appears at control line (C) within 15 minutes of performing the test, then the test has not worked. You should test again, using a new test, taking care to follow the instruction. At the same time, contact our email: bio@micropf.com immediately.

Q10. Can any medication or medical conditions affect the results?

We have done research on the effects of the drug, see the **INDEX OF CHARACTERISTICS**. The results showed that the drugs listed in Section 5 had no effect on the test results. If you are taking medicines other than those listed, ask your doctor for advice.

Q11. What are the possible risks of this test?

Possible Risks:

- Discomfort during the sampling
- Incorrect test results (see Interpreting Results and Limitations Sections).

Q12. What should I do if there is blood on the saliva ?

Please check your physical condition. If the bleeding is serious, you should contact the doctor immediately. Blood will not affect the test results.

INDEX OF CHARACTERISTICS

1. Positive reference coincidence rate: the positive reference coincidence rate of the enterprise should be 100%.
2. Negative reference product conformity rate: the negative reference product conformity rate of the enterprise should be 100%
3. Limit of detection (LoD):The LoD is determined using limiting dilutions of inactivated SARS- CoV-2 in two separate methods. The inactivated virus is spiked into the treatment buffer processed with a negative saliva sample to have a concentration of TCID₅₀ of 5.6 x 10⁵/mL. Each sample is serially 10-fold diluted and by testing in triplicate, a tentative LoD showing 100% (3/3) positive rate is determined for each. For confirmation LoD study, 4 concentrations below the lowest concentration of the pre-test are tested in 20 replicates and a concentration showing over 95% (19/20) are positive, determined as the LoD of the fluorecare® SARS-CoV-2 Antigen Test Kit. This was: 49 TCID₅₀/mL.
4. Cross-reactivity: Virus/bacteria listed below are confirmed not to have cross-reactivity with SARS-CoV-2 Antigen Test Kit.

Human Coronavirus (OC43) 3.8 x 10⁵PFU/ml ;Human Coronavirus (229E) 2.3 x 10⁴PFU/ml ;Human Coronavirus MERS (Florida/USA-2_Saudi Arabia_2014) 1.05 x 10⁵PFU/ml ;Human Coronavirus (NL63) 2.8 x 10⁴PFU/ml ;Human Coronavirus (HKU1) (N-protein) 45µg/ml ;Adenovirus Type 01 (Species C) 8.34x 10⁴PFU/ml ;Adenovirus Type 02 (Species C) 1.05 x 10⁵PFU/ml ;Adenovirus Type 11 (Species B) 1.02 x 10⁷PFU/ml ;Enterovirus Type 68 (2014 Isolate) 1.05 x 10⁵PFU/ml ;Human Metapneumovirus(16 Type A1) 3.80 x 10⁵PFU/ml ;Human Metapneumovirus (3 Type B1 strain Peru 2_2002) 1,41 x 10⁴PFU/ml ;Parainfluenza Virus (Type 1) 1.26 x 10⁵PFU/ml ;Parainfluenza Virus (Type 2) 1.26 x 10⁵PFU/ml ;Parainfluenza Virus (Type 3) 3.39 x 10⁶PFU/ml ;Parainfluenza Virus (Type 4B) 3.80 x 10⁵PFU/ml ;Respiratory Syncytial Virus Type A (Isolate: 2006) 7.35 x 10⁵PFU/ml(RSV excepted) ;Rhinovirus (Type 1A) 1.05 x 10⁶PFU/ml ;Influenza Type A, H3N2

(HK/8/68) 1.51 x 10⁴PFU/ml(Influenza A excepted) ;Influenza Type A, H1N1 (Brisbane/59/07) 4.57 x 10⁵PFU/ml(Influenza A excepted) ;Influenza Type A, H1N1pdm (Canada/6294/09) 1.26 x 10⁵PFU/ml(Influenza A excepted) ;Influenza Type B (Texas/6/11) 2.26 x 10⁵PFU/ml(Influenza B excepted) ;Influenza Type B (Alabama/2/17) 3.16 x 10⁵PFU/ml(Influenza B excepted) ;Staphylococcus aureus (Protein A) DSM 21705 (E. Domann) 3.62 x 10⁹CFU/ml ;Staphylococcus aureus (Protein A) DSM 21979 (E. Domann, Univ.) 7.64x 10⁹CFU/ml ;Staphylococcus aureus (Protein A) DSM 46320 (E. Domann) 4.58x 10⁹CFU/ml ;Staphylococcus epidermidis DSM 1798 (PCI 1200) 4.90 x 10⁹CFU/ml ;Staphylococcus epidermidis DSM 20044 (Fussel) 5.10 x 10⁹CFU/ml ;Bordetella pertussis DSM 4923 (Walker) 2.71 x 10⁹CFU/ml ;Bordetella pertussis DSM 4926 (Sato and Arai) 2.02 x 10⁹CFU/ml ;Bordetella pertussis DSM 5571 8.07 x 10⁹CFU/ml ;Legionalle pneumophila DSM 7513 (Philadelphia-1) 4.50 x 10⁹CFU/ml ;Legionalle pneumophila DSM 7514 (Los Angeles-1) 1.17 x 10¹⁰CFU/ml ;Streptococcus pyogenes DSM 20565 (SF130, T1) 1.37 x 10⁸CFU/ml ;Streptococcus pyogenes DSM 2071 (S. Koshimura, Sv) 9.30 x 10⁷CFU/ml ;Haemophilus influenzae DSM 24049 (TD-4) 7.77 x 10⁸CFU/ml ;Haemophilus influenzae DSM 4690 (Maryland) 1.41 x 10⁷CFU/ml ;Haemophilus influenzae DSM 23393 (Pittman 576) 1.23 x 10⁸CFU/ml ;Mycobacterium tuberculosis DSM 43990 (BCGT, tice) 4.69 x 10⁸CFU/ml ;Streptococcus pneumoniae (Protein G) DSM 20566 (SV1) 4.05 x 10⁸CFU/ml ;Streptococcus pneumoniae (Protein G) DSM 11967 (Jorgensen262) 3.80 x 10⁷CFU/ml ;Streptococcus pneumoniae (Protein G) DSM 25971 (Gyeonggi) 2.70 x 10⁸CFU/ml ;Mycoplasma pneumoniae DSM 23978 (Eaton Agent, FH) >10⁵ cells/ml ;Mycoplasma pneumoniae DSM 23979 (M129-B7) >10⁵ cells/ml ;Candida albicans DSM 1386 (NIH 3147) 6.53 x 10⁵CFU/ml ;Candida albicans DSM 1665 (132) 2.39 x 10⁸CFU/ml ;Candida albicans DSM 5817 (806M) 2.55 x 10⁸CFU/ml ;Pseudomonas aeruginosa DSM 1117 (Boston 41501) 1.31 x 10⁹CFU/ml ;Pseudomonas aeruginosa DSM 3227 (Schutze) 3.93 x 10⁹CFU/ml ;Streptococcus salivarius DSM 20560 (275) 5.44 x 10⁸CFU/ml ;Streptococcus salivarius DSM 20067 (21367) 5.09 x 10⁸CFU/ml ;

5. Interference

Substances listed below are confirmed not to have interference response with SARS-CoV-2 Antigen Test Kit.

Benzocaine (150 mg/dL), Blood (human) (5%), Mucin(5 mg/mL), Naso GEL (NeilMed) (5%), CVS Nasal Drops (phenylephrine) (15%), Afrin (Oxymetazoline)(15%),CVS Nasal Spray(Cromolyn(15%),Zicam Cold Remedy (5%), Homeopathic (Alkalol) (10%), Sore Throat Phenol Spray (15%),Tobramycin(3.3mg/dL),Mupirocin(0.15mg/dL),Fluticasone (0.000126 mg/dL),Tamiflu (Oseltamivir phosphate) (500mg/dL), Budenoside (0.00063 mg/dL), Biotin (0.35mg/dL),Methanol(150mg/dL),Acetylsalicylic Acid(3mg/dL), Diphenhydramine (0.0774mg/dL),Dextromethorphan (0.00156mg/dL), Dexamethasone (1.2 mg/dL), Mucinex(5%).

6.Hook

Virus strains	Limit value
SARS-CoV-2	5.6 x 10 ⁵ TCID ₅₀ /mL

7. Clinical Accuracy

Test method: RT-PCR positive participants within 7 days after symptoms were tested and the positive coincidence rate of the kit was compared. Number of test samples:818 samples, including 368 RT-PCR positive samples; 335 SARS-CoV-2 Antigen Test Kit positive results;33 SARS-CoV-2 Antigen Test Kit negative results;450 RT-PCR negative samples; 450 SARS-CoV-2 Antigen Test Kit negative results.

Samples: In total, 818 participants are included for the unit and all test samples included are tested.

Compared with self-test and RT-PCR Test, 32 positive RT-PCR samples; 30 SARS-CoV-2 Antigen Test Kit results were positive, 2 SARS-CoV-2 Antigen Test Kit negative results; 87 RT-PCR negative samples, 87 SARS-CoV-2 Antigen Test Kit negative results

6.1.Statistics on test results and those of the product tested are as follows:

Method	RT-PCR		Total Results
	Results	Positive	
SARS-CoV-2 Antigen Test Kit (Colloidal Gold Chromatographic Immunoassay)(saliva)	Positive	335	335
	Negative	33	450
	Total Results	368	818

Method		RT-PCR		Total Results
by Lay person	Results	Positive	Negative	
	Positive	30	0	30
	Negative	2	87	89
Total Results		32	87	119

6.2 PCR Cycle Threshold (CT) analysis of PPA at each stage

Cycle Threshold (CT)	# of RT-PCR positive	fluorecare® SARS-CoV-2 Antigen Test (Colloidal Gold Chromatographic Immunoassay)		
		# of positive results	PPA	NPA
<25	105	104	99.05%	100%
<30	217	214	98.62%	
<35	297	292	98.32%	
<38	368	335	91.03%	

Clinical sensitivity at Ct<25 = 104/105 = 99.05%

Positive correct rate(Clinical sensitivity) at Ct<38=335/368=91.03%

(95%CI:87.64%-93.75%)

Negative correct rate(Clinical specificity) = 450/450 = 100% (95%CI:99.18%~100.00%)

The result of fluorecare® SARS-CoV-2 Antigen Test (Colloidal Gold Chromatographic Immunoassay): clinical sensitivity was 98.62% (Ct <30) and the clinical specificity was 100%

7. Repeatability: The repeatability reference products of the enterprise were tested, repeated for 10 times, and the positive coincidence rate is 100%.

8. The fluorecare® SARS-CoV-2 Antigen Test Kit (Colloidal Gold Chromatographic Immunoassay) is test for SARS-CoV-2 nucleocapsid protein, the mutants of SARS-CoV-2 Alpha, Beta, Gamma, Delta and Omicron can be identified by the fluorecare® SARS-CoV-2 Antigen Test Kit.

WARNING AND PRECAUTION

1. Read the Instruction for use completely before using the product. Follow the instructions carefully. Failure to do so may result in an inaccurate result.
2. The kit is only used for in vitro diagnosis; it cannot be used repeatedly. Do not swallow.
3. Avoid getting the buffer solution into the eyes or skins.
4. Keep out of reach children.
5. The test kit is for single use only, do not reuse any components of the test kit.
6. Do not use this test beyond the expiration date printed on the outer package. Always check expiry date prior to testing.
7. Do not touch the reaction area of the test cassette.
8. Do not use the kit if the pouch is punctured or not well sealed.
9. DISPOSAL: All specimens and the used-kit has the infectious risk. The process of disposing the diagnostic kit must follow the local, state and federal infectious disposal laws/regulations.
10. During the time of interpretation, no matter the shade of the color band, it can be found to be positive as long as two lines appear on the quality control area and the detection area, respectively.
11. Please ensure that an appropriate amount of sample is used for testing, too much or too little of sample amount will cause the result deviation.
12. The final result should be read in 15 minutes. Please do not read the result after 20 minutes.
13. Various components of different batch of reagents cannot be used interchangeably in order to avoid wrong results
- 15.If the terminal on the top of the cap is not unscrewed, and if the blue stopper inside the sample processing tube is not removed, it will not be possible to drip liquid!
16. Only 2 drops of the treated sample solution can be added ! Adding too much or too little of the treated sample solution may result in invalid test results.

INTERPRETATION OF ICONS

	Do not re-use		Temperature limit
	In vitro diagnostic medical device		Consult instructions for use
	Contains sufficient for <n> tests		Authorized representative in the European Community
	Keep away from sunlight		Caution
	Manufacturer		CE marking
	Catalogue number		Batch code
	Date of manufacture		Use-by date
	Self-test		Do not use if package is damaged

GENERAL INFORMATION

Shenzhen Micropofit Biotech Co., Ltd.

Rm. 405, 406, Zone B /4F, Rm. 205, 206-1, 207, West Side of Zone B/ 2F, Haoui Building, No. 8 Langshan 2nd Road, Songpingshan, Songpingshan Community, Xili Street, Nanshan District, Shenzhen, P.R. China

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